

MRI micro- and macro- structural biomarkers linked to nonpathological brain aging in macaques

L. Dieckhaus¹, K. McDermott^{2,4}, A. Murlikrishnan³, E. M. Foiryolo³, D. T Gray⁴,
E. Hutchinson¹, C. A Barnes^{2,4}



¹Department of Biomedical Engineering, University of Arizona ²Departments of Psychology, Neurology
And Neuroscience, University of Arizona, ³KEYS Res. Internship, BIO5 Inst., ⁴Evelyn F. McKnight Brain
Institute

Society for Neuroscience 2021



THE UNIVERSITY OF ARIZONA
COLLEGE OF ENGINEERING


Biomedical Engineering

Introduction

Using a robust Magnetic Resonance Imaging (MRI) microscopy battery composed of

- Diffusion MRI (dMRI)
- Relaxometry based MRI (rMRI)

Goal: To quantify changes in brain structure related to normal aging.




- Cohort Size: 12 female non-human primates (bonnet macaques) brains provided by the Barnes lab.
 - Ages ranging from 10 to 25 years old (30 to 75 human equivalent years).
- Ex-vivo MRI battery is developed to examine brain specimens at high resolution (**200-600 micron scale**) for quantitative biomarkers that indicate age-related brain dysfunction.

Background

Multiple Spin Echo (MSE) – rMRI Technique:


Myelin Water Fraction (MWF) – Quantifies water trapped between myelin bilayers (a marker for myelin content – quantifies myelin found in each voxel).



Selective Inversion Recovery (SIR)

– rMRI Technique:


Bound Pool Fraction (BPF) – Measures myelin content and is sensitive to molecular size.



Diffusion Tensor Imaging (dMRI):


Fractional Anisotropy (FA) – Reports that there is a preferred direction of water diffusion in each voxel.

Trace (TR) – Reports the total water diffusion regardless of direction.




Summary of Methods

Sample Preparation



Acquisition: Ex- Vivo Battery

MR Battery (~50 hours)



Map Generation


Different pipelines avenues depending on scan type.



Summary of Methods


Image Registration

Using Advanced Normalization Tools (ANTs) to register native space images to template space.





Voxel-wise Analysis

Using fsl randomise tool threshold free cluster enhancement and assessing negative and positive correlations with age.




ROI Analysis

Masked high resolution anatomical (HRA) and conducted a Pearson correlation test across cohort.




Results – Voxel-wise Analysis


Fractional Anisotropy




Trace



Bound Pool Fraction



Myelin Water Fraction




Threshold Free Cluster
Enhancement corrected p value
5000 permutations (conservative)




Results – ROI Analysis


- Found strong correlations when assessing interactions between different MR metrics to one another.




Hippocampus: FA vs. MWF




Hippocampus: MWF vs. BPF




Hippocampus: FA vs. BPF





White Matter: FA vs. MWF



White Matter: FA vs. BPF



White Matter: MWF vs. BPF



Conclusions

In Voxel-wise and ROI analysis with respect to aging:

- No significant positive or negative correlations were found between our MR metrics and age using voxel-wise analysis.
- No significant positive or negative correlations were found in the hippocampus or white matter with respect to age using ROI analysis.
- **Both analyses produced consistent results with regards to aging.**


In ROI analysis with respect to MR metric correlation:

- Strong positive correlation with FA and MWF measures in white matter were observed.
- Strong positive correlation with BPF and MWF measures in the hippocampus.
- This indicates that these 3 measures, depending on location in the brain, contain similar information.
 - The combination of MR metrics may be better at identifying anatomical information than a single metric on its own.

Future Directions


- Quantitative Susceptibility Maps (QSM) which is sensitive to iron and calcium accumulation.
 - Voxel-wise and ROI analysis with behavioral data and age.
- Use principal component analysis to perform dimensionality reduction to test if the combination multiple MRI metrics (features) can characterize normal aging.
 - ROI analysis of subfields of the hippocampus.

Acquired Echoes (Magnitude & Phase)




N = 1 - 4
TE = 4.5 - 28.1 ms


QSM



Behavior



MR Metrics



References

REMMI-matlab (<https://github.com/remmi-toolbox/remmi-matlab>)

TORTOISE (<https://tortoise.nibib.nih.gov/>)

ANTs (<https://github.com/ANTsX/ANTs>)

Lowekamp BC, Chen DT, Ibáñez L and Blezek D (2013) The Design of SimpleITK. *Front. Neuroinform.* 7:45. doi: 10.3389/fninf.2013.00045

Z. Yaniv, B. C. Lowekamp, H. J. Johnson, R. Beare, "SimpleITK Image-Analysis Notebooks: a Collaborative Environment for Education and Reproducible Research", *J Digit Imaging.*, <https://doi.org/10.1007/s10278-017-0037-8>, 2017.

Winkler AM, Ridgway GR, Webster MA, Smith SM, Nichols TE. Permutation inference for the general linear model. *NeuroImage*, 2014;92:381-397.

Acknowledgements

The authors would like to thank the Machining & Welding Center at Gould Simpson, in particular Larry Acedo, Biomedical Device Prototyping Services, in particular, Scott DeRigne, 3D Printing Facilities, in particular Charles Ingram, 7T and Translational Imaging Research, the High Performance Computing (HPC) Core, and Biorender which was utilized for creating the animated images. Also, the authors would like to thank Avinash Murlikrishnan and Eric Moses Foiryolo for the outstanding work on the masks that drove the ROI preliminary analysis. The main author would like to thank Christa Sonderer for helping with obtaining preliminary data from the QSM processing pipeline. The authors would also like to thank Arne Ekstrom for providing a human MRI hippocampal subfield segmentation method that we are currently adhering to the bonnet macaque cohort.